



# ASPIRATION - IRRIGATION SYSTEMS







## ASPIRATION HANDPIECE




LENGTH	Ga	CODE
105 mm	21	TMN 213
105 mm	22	TMN 103
105 mm	23	TMN 113

## IRRIGATION HANDPIECE



LENGTH	Ga	CODE
105 mm	21	TMN 214
105 mm	22	TMN 104
105 mm	23	TMN 114


## ASPIRATION CANNULA, 22 Ga




TMN 107




## IRRIGATION CANNULA, 22 Ga



TMN 108



ASPIRATION HANDLE  
TMN 105





IRRIGATION HANDLE  
TMN 106

MONOZOOM

SEPARABLE


# COAXIAL A/I SYSTEM



COAXIAL I/A SYSTEM  
USED WITH VARIOUS SILICONE  
SLEEVES FROM 1.8 MM INCISION

ASPIRATION CANNULA  
FOR AIS, STRAIGHT, 22 Ga

TMN 201




ASPIRATION CANNULA  
FOR AIS, CURVED, 22 Ga

TMN 202



ASPIRATION CANNULA  
FOR AIS, STRAIGHT, 22 Ga  
WITH BALL 0,8 mm DIA.

TMN 203



ASPIRATION CANNULA  
FOR AIS, CURVED, 22 Ga  
WITH BALL 0,8 mm DIA.

TMN 204



E  
R  
G  
O  
N  
O  
M  
I  
C  
  
H  
A  
N  
D  
L  
E





Technicheskaya str. 120 A  
420054, P.O. Box 278, Kazan, Russia

Phone: +7 917 2610893  
Phone / Fax: +7 843 277 07 78, 260 17 58  
+7 843 558 07 78  
[www.titanmedical.ru](http://www.titanmedical.ru)  
E-mail: [sales@titanmedical.ru](mailto:sales@titanmedical.ru)